

Adox

Oral Care Series

mouth+aid



Medical Care Labo



Oral health

Adox

mouth+aid

For oral care while sleeping.

Adox.mouth+aid uses bacteriostatic materials.
This material is patented by us and the effect
The effect is long-lasting.

Effectiveness of the mouthpiece.

The following have been proven simply by wearing it during sleep.

1 Oral periodontal disease, bacteria reduced

According to a paper published by Kanagawa Dental University,
it can be worn while sleeping.

The oral cavity is kept in a healthy state by inhibiting the growth of
periodontal disease bacteria in the mouth.

2 Decrease in methylcaptan.

It reduces methylcaptan, which can also cause bad breath.

Recommended for people who are concerned about their own bad breath.

3 Antibacterial effect on mouthpieces.

This mouthpiece has an added antimicrobial ingredient that The antimicrobial
effect lasts longer.

Development/Raw materials. : Adox Inc.

Manufactured/distributed :Medical Care Labs Inc.

URL : adox.pro/pd

(Head office) 2-19-1 Nishiki, Naka-ku, Nagoya
460-0003, Nagoya Konoike Building 1104
(Research Institute) 1-23-2, Mizuhodori, Mizuho-ku,
Nagoya 467-0806, King Heights Mizuho 2F
TEL (052) 841-3677 FAX (052) 841-9889



Adox

Oral Care Series

mouth+aid


MouthPiece

For oral care while sleeping....

Adox.mouth+aid uses a bacteriostatic material. This material is patented by our company and has a lasting effect.


Easy to mould to your own tooth mould.

1




Place the mouthpiece in hot water at approx. 75 °C for approx. 30 seconds. Remove with chopsticks or similar to avoid scalding.

2




In front of a mirror, bite down firmly in front of the mirror in a position where the indentation of the mouthpiece is centred on the upper teeth.
※Beware of burns in the ministry.

3



Place your fingers on the outside of the upper teeth and press them against the teeth and gums so that they fit together.

4



Mouthpiece. Remove it from your mouth. Cool it with water. The process is complete.

- Product name : Adox mouth+aid ■ Material: EVA resin body Case: PP resin
- Product size Width 5.3cm x Length 6.3cm x Height 1.7cm
- Contents: Mouthpiece x 1, Japanese instructions, special case.
- ※ Specifications and product design are subject to change without notice.

◆ Storage precautions

- After use, wash with AdoxPG2500 antibacterial cleaning agent or denture cleaner and dry well before Store in the case.
- Do not store in extremely high or low temperatures or in direct sunlight.
- Keep out of reach of children.

〈Intended use.〉

Do not use the product for any purpose other than its intended use.

Manufacturer's suggested retail price [incl. tax] 19,800 yen

◆ Read before use.

〈Precautions for use.〉

The mouthpiece overlaps the upper teeth and swells around the mouth, causing discomfort in the mouth. Please start with a few hours during the day, and then start using the mouthpiece during sleep until you feel no more discomfort.

〈Do not use if you have any of the following.〉

- People with breathing difficulties due to nasal disorders.
- People who are undergoing dental treatment or have oral abnormalities.
- People with many dentures.
- If you experience pain in the temporomandibular joint when using the device, please stop using it.
- When using the device on children, please consult a dentist or other specialist before use.

Development/raw materials : ADOX Co.
Manufactured/distributed: Medical Care Labo Co.
(Head office) 2-19-1 Nishiki, Naka-ku, Nagoya 460-0003, Nagoya Konoike Building 1104
(Laboratory) 1-23-23, Mizuho-dori, Mizuho-ku, Nagoya 467-0806, King Heights Mizuho 2F,
TEL (052) 841-3677 FAX (052) 841-9889
homepage <http://mcl.co.jp>

Sole distributor [Dealer/dealer]